

RANGE

After more than 3 decades of supremacy in the Indian market, the Kiato range of surgical blades, disposable scalpels and stitch cutters made their global presence felt in 2000. Today, over 60 countries bear testimony to the excellent quality and superior performance of our products, both under the Kiato brand name as well as in private label. While attempting to realise our dream of becoming the largest and finest producer of surgical blades in the world, we are also providing excellent value for money.

By creating innovative new products and improving the quality of our existing portfolio, we consistently strive for excellence. Towards this end, we have upgraded our ISO 9002 registration to ISO 9001 in 2001. The overall product quality, with special emphasis on sharpness, has been improved by over 30% and independent tests have confirmed that the sharpness of our blades is better than most leading blades. Podiatry blades have also been added to our range.

Medica 2002 will witness the launch of a new safety scalpel offering better features and good rigidity and will reinforce our attempts to be at the forefront of developments within the industry.



Mini Scalpel
ABS plastic handle fitted with a stainless steel blade. Protected with a guard. Ideal for kits.



Disposable Scalpel
Sturdy ABS plastic handle with removable plastic guard provides protection. Sizes 10, 11, 12, 15, 20, 21, 22, 23 & 24.



Safety Scalpel

Sterile Disposable safety scalpel with retractable sheath. The right choice to help you meet OSHA guidelines. Prevents sharp related injuries. Available in sizes 10, 11 & 15.



Metal BP Handle
Individually packed. Fit our blades perfectly. Made of stainless steel.



Kiato



CE
0 4 3 4

RANGE

S
e
d
e
s

Podiatry Blades

Individually packed in foil. Packed 12 blades per box. Sterile.



Stitch Cutters

Available in 2 styles. Suitable for No. 3 Fitment



Gouge Blades

Very sharp edge made of carbon steel. Individually packed in foil.



Surgical Blades

Offered in carbon steel as well as stainless steel.



Fitment No:3

Fitment No:4

Plot No.19 Sector 6 I.M.T. Manesar
Gurgaon Haryana 122 050 India
T: +91 124 2291406/7 F: +91 124 2291206
Email: info@adityadispomed.com
www.adityadispomed.com

aditya
DisPomed
PRODUCTS
PRIVATE
LTD