



*The smart choice for sharp blades*



## Premium Sterile Surgical Blade

- Carbon Steel Blades made of Swiss Steel grade 125Cr1.
- Stainless Steel Blades made of premium Swedish Steel.
- Individually packed in easy peel color coded foil and boxes.
- Fits all makes of handles. Conforms to ISO 7740.
- Sterilized by Gamma Radiation to a minimum dosage of 2.5 MRads.
- Blades conform to BS 2982 and Medical Device Directive 93/42/EEC.



100 blades per box, 10 boxes per case, 10 cases per master shipping carton.



Size	10	10A	11	11/11P	12	12D	15	15/15T	15C
Carbon Steel	18010	18009	18011	18031	18012	18025	18035	18015	18026
Stainless Steel	18110	18109	18111	18131	18112	18125	18135	18115	18126



Size	18	20	21	22	22A	23	24	25	36
Carbon Steel	18027	18020	18021	18022	18030	18023	18024	18028	18029
Stainless Steel	-	18120	18121	18122	-	18123	18124	18128	-

Blade size not to scale.

All products are CE marked as per the European Medical Device Directive 93/42/EEC.



EDITION 02/May 2010



0 4 3 4

aditya  
DISPOMED



Plot No. 19 Sector 6 IMT Manesar  
Gurgaon Haryana I22050 India

Tel. +91 I24 22 91 406 / 7

Fax +91 I24 22 91 206

info@adityadispomed.com

www.adityadispomed.com

A subsidiary of Kehr Surgical Pvt. Ltd.