

Kiato⁺ PLUS



The smart choice for sharp blades

- Carbon Steel Blades made of Swiss steel grade 125Cr1.
- Stainless steel blades made of premium Swedish Steel.
- Individually packed in easy peel color coded foil and boxes.
- Fits all makes of handles. Conforms to ISO 7740.
- Sterilized by Gamma Radiation to a minimum dosage of 2.5 MRads.
- Blades conform to BS 2982 and Medical Device Directive 93/42/EEC.



STEEL SPECIFICATION

Carbon Steel	
Constituent	Percentage
Carbon	1.20 to 1.30
Silicon	0.10 to 0.35
Manganese	0.10 to 0.45
Chromium	0.10 to 0.40
Sulphur	0.025 Max.
Phosphorus	0.035 Max.
Easy Peel with protective German VCI paper	

Stainless Steel	
Constituent	Percentage
Carbon	0.6 to 0.7
Manganese	1.0 max.
Chromium	12.0 to 13.5
Silicon	0.5 max.
Phosphorus	0.03 max.
Sulphur	0.025 max.
Nickel	0.5 max.

100 blades per box, 10 box per case, 10 case per master shipping carton.

The range of specialty blades finds application in Dentistry, Histopathology and even for Dermaplaning besides general surgery. Our blades are made with Passion using our experience of 50+ years. The combination of premium inputs and equipment from Japan, Germany, Switzerland, France and UK is your assurance for reliability.

- German Anti Corrosive Paper for protection
- Laser controlled edge for sharpness & Durability



100 blades per box, 10 box per case, 10 case per master shipping carton.

FITMENT 3



Size	10	11	11/11P	12	15	15/15T
Carbon Steel	18010	18011	18031	18012	18035	18015
Stainless Steel	18110	18111	18131	18112	18135	18115

Blade size not to scale.

FITMENT 4



Size	19	20	21	22	23	24
Carbon Steel	18019	18020	18021	18022	18023	18024
Stainless Steel	18119	18120	18121	18122	18123	18124

Special Sizes

FITMENT 3



Size	10A	10R	12D	13	14	15C
Carbon Steel	18009	18040	18025	18013	18014	18026
Stainless Steel	18109	18140	18125	18113	18114	18126

Blade size not to scale.

FITMENT 4



Size	18	22A	25	26	36
Carbon Steel	18027	18030	18028	18033	18029
Stainless Steel	18127	-	18128	18133	18129